visionrt

MapRT DICOM Conformance Statement

For MapRT v1.0, v1.1, v1.2, v1.2.1

1016-0591-EN October 2024

Document Version: 1.2



Vision RT Ltd. Dove House Arcadia Avenue London N3 2JU UK

EN-US



Information in this document is subject to change without notice. Stenographical and clerical errors on the face of this document are subject to correction. By accepting and using this document, the Customer agrees that Vision RT's liability for errors contained in this document or loss or damage incurred by the Customer in connection with the furnishing or use of this document shall be limited to and in accordance with the warranty given to the product(s) supplied, as set out in Vision RT's standard terms and conditions of sale, which apply to and govern this document. Vision RT shall not be responsible for any other loss, liability, or damages whatsoever. This document may contain information confidential and/or proprietary to Vision RT including, without limitation, technical or commercial know-how, specifications, inventions, processes, or initiatives which are of a confidential nature. This document and its contents must therefore be kept in strict confidence by the Customer in accordance with Vision RT's standard terms and conditions of sale and may only be used for the purpose for which it is intended. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying and recording, for any purpose without the express written permission of Vision RT Ltd. Nothing in this document shall constitute the granting (expressly or impliedly) of any right or license to the Customer to any intellectual property contained herein.

Trademarks All other products or services mentioned in this document are identified by the trademarks or service marks of their respective companies or organizations. Vision RT Ltd disclaims any responsibility for specifying which marks are owned by which companies or organizations.

Copyright © 2024 Vision RT Ltd.

All rights reserved.

1 OVERVIEW

MapRT is a treatment delivery simulation system used to aid the treatment planning process and alert the user to potential collisions with equipment or patient. To support this, it acts as a provider of service supporting transfer of RT Plans.

This document explains the assumptions, limitations, and supported format.

SOP Classes		User of	Provider	
Name	UID	(SCU)	(SCP)	
Transfer				
RT Plan Storage	1.2.840.10008.5.1.4.1.1.481.5	No	Yes	

Figure 1 – Network Services Supported by MapRT



2 TABLE OF CONTENTS

1	OVERV	IEW	3
2	TABLE	OF CONTENTS	4
3	INTRO	DUCTION	6
	3.1 Re	vision history	6
	3.2 Au	dience	6
	3.3 Re	marks	6
	3.4 De	finitions, Terms, and Abbreviations	6
4	NETWO	DRKING	9
	4.1 Im	plementation Model	9
	4.1.1	Application Data Flow	9
	4.1.2	Functional Definition	9
	4.1.3	Sequencing of Real-World Activities	9
	4.2 Spe	ecification	10
	4.2.1	SOP Classes	
	4.2.2	Association Policies	
	4.2.3	Association Initiation Policy	
	4.2.4	Association Acceptance Policy	
	4.3 Ne	twork Interfaces	11
	4.3.1	Physical Network Interface	
	4.3.2	Additional Protocols	
	4.3.3	IPv4 and IPv6 Support	
	4.4 Co	nfiguration	11
	4.4.1	AE Title/Presentation Address Mapping	
	4.4.2	Parameters	11
5	MEDIA	INTERCHANGE	12
6	TRANS	FORMATION OF DICOM TO CDA	13
7	SUPPO	RT OF CHARACTER SETS	14
, ,	SECURI		15
ð			15 15
		control lovel Socurity	15
	0.2 AS: 9.2 An	nlication Level Security	15
_	0.5 AP		
9	ANNEX	ES	16
	9.1 101	O Contents	16
	9.1.1	Created SOP Instance(s)	
	9.1.2	Usage of Attributes from Received IODs	16
	9.1.3	Attribute Mapping	
	9.1.4	Coerced/Modified Fields	
	9.2 Da	ta Dictionary of Private Attributes	
	9.3 Co	aed Terminology and Templates	

visi∳∩**rt**

9.4	Grayscale Image Consistency	19
9.5	Standard Extended/Specialized/Private SOP Classes	19
9.6	Private Transfer Syntaxes	19

3 INTRODUCTION

3.1 Revision history

Date	Rev	Author	Comments
22 Jun 2023	1.0	Andrzej Wawrzyńczyk	Initial version for MapRT 1.0 using the NEMA template at: <u>https://dicom.nema.org/medical/dicom/2022d/output/pdf/part02.pdf</u> (Chapter A)
16 Feb 2024	1.1	Michelle Kennedy	This document is applicable to MapRT v1.1 as well as v1.0
17 Oct 2024	1.2	Andrzej Wawrzyńczyk	This document is applicable to MapRT v1.2 and v1.2.1.

3.2 Audience

The target audience are hospital IT support staff, and service and support engineers of other radiotherapy product vendors, who need to configure DICOM interoperability between MapRT and their products.

It is assumed that the reader is familiar with the DICOM standard.

The document is structured based on the template definition in Part 2 of the DICOM standard (edition 2022d).

3.3 Remarks

It is essential that this DICOM Conformance Statement is thoroughly reviewed when using MapRT with other DICOM capable systems.

3.4 Definitions, Terms, and Abbreviations

Informal definitions are provided for the following terms used in this Conformance Statement. The DICOM Standard is the authoritative source for formal definitions of these terms, and is available for download at https://www.dicomstandard.org/current

Abbreviation/Term	Explanation
Abstract Syntax	The information agreed to be exchanged between applications, generally equivalent to a Service/Object Pair (SOP) Class. Examples: Verification SOP Class, Modality Worklist Information Model Find SOP Class, Computed Radiography Image Storage SOP Class.
AE	Application Entity
	An end point of a DICOM information exchange, including the DICOM network or media interface software. This may be software that sends or receives DICOM information objects or messages. A single device may have multiple AEs.
AET	Application Entity Title
	The externally known name of an <i>Application Entity</i> , used to identify a DICOM application to other DICOM applications on the network.

Abbreviation/Term	Explanation
Attribute	A data element identified by a Tag in an object definition. This may itself be a complex data structure (Sequence) composed of lower-level data elements.
	Examples: Patient ID (0010,0020), Accession Number (0008,0050), Photometric Interpretation (0028,0004), Procedure Code Sequence (0008,1032).
DICOM	Digital Imaging and Communications in Medicine
	A standard for the management and transmission of medical images and related data used in many healthcare facilities.
IOD	Information Object Definition
	The specified set of Attributes that comprise a type of data object; does not represent a specific instance of the data object, but rather a class of similar data objects that have the same properties.
	The Attributes may be specified as Mandatory (Type 1), Required but possibly unknown (Type 2), or Optional (Type 3), and there may be conditions associated with the use of an Attribute (Types 1C and 2C). Examples: MR Image IOD, CT Image IOD, Print Job IOD.
Module	A set of Attributes that are logically related to each other within an IOD.
	Example: Patient Module includes Patient Name, Patient ID, Patient Birth Date, and Patient Sex.
NEMA	National Electrical Manufacturers Association
	The standards body that publishes the DICOM standard; <u>https://www.nema.org/</u>
PACS	Picture Archiving and Communication System
	A medical imaging technology used primarily in healthcare organizations to securely store and digitally transmit electronic images and clinically relevant reports.
RT	Radiotherapy / Radiation Therapy
	A treatment that uses high-energy x-ray or other particles to destroy cancer cells.
SCP	Service Class Provider
	The role of an Application Entity that provides a DICOM network service; typically, a server that performs operations requested by another Application Entity. Examples: Picture Archiving and Communication System (image storage SCP, and image query/retrieve SCP), Radiology Information System (modality worklist SCP).

Abbreviation/Term	Explanation	
SCU	Service Class User	
	The role of an Application Entity that uses a DICOM network service; typically, a client. Examples: imaging modality (image storage SCU, and modality worklist SCU), imaging workstation (image query/retrieve SCU)	
Sequence	A structure of repeated data within a DICOM file. Members of a sequence are identified using the symbol '>' before their Attribute or Module names.	
SOP	Service Object Pair	
	Specifies the network or media transfer (service) of a particular type of data; the fundamental unit of DICOM interoperability specification. A valid DICOM file shall have one (and only one) SOP class UID assigned to it.	
	Examples: Ultrasound Image Storage Service, Basic Grayscale Print Management.	
Тад	Data Element Identifier	
	A 32-bit identifier, represented as a pair of four-digit hexadecimal numbers describing the "group" and the "element". If the "group" number is odd, the tag is for a private (manufacturer-specific) data element.	
	Examples: (0010,0020) [Patient ID], (07FE,0010) [Pixel Data], (0019,0210) [private data element].	
TPS	Treatment Planning System	
	A computer system aiding the radiation treatment planning. E.g., Varian Eclipse, Philips Pinnacle, RaySearch RayStation, Elekta Monaco.	
UID	Unique Identifier	
	A globally unique "dotted decimal" string that identifies a specific object instance or a class of objects; an ISO-8824 Object Identifier.	
	Examples: Study Instance UID, SOP Class UID such as 1.2.826.0.1.3680043.2.120.20000624	
VR	Value Representation	
	The format type of an individual DICOM data element, such as text, an integer, a person's name, or a code. DICOM information objects can be transmitted with either explicit identification of the type of each data element (Explicit VR), or without explicit identification (Implicit VR); with Implicit VR, the receiving application must use a DICOM data dictionary to look up the format of each data element.	

4 NETWORKING

4.1 Implementation Model

4.1.1 Application Data Flow



4.1.2 Functional Definition

The MapRT DICOM Listener application entity is implemented as a Windows service. It is started at server start-up and accepts storage and verification requests from remote DICOM nodes. It acts as a service class provider (SCP) for C-ECHO and C-STORE requests. It only accepts an incoming connection and never initiates a new connection.

4.1.3 Sequencing of Real-World Activities

Not applicable.

4.2 Specification

4.2.1 SOP Classes

The MapRT DICOM Listener provides Standard Conformance to the following class(es):

S	SOP Classes		Provider	
Name	UID	(SCU)	(SCP)	
RT Plan Storage	1.2.840.10008.5.1.4.1.1.481.5	No	Yes	

4.2.2 Association Policies

4.2.2.1 General

The DICOM standard application context shall be used:

Name	UID
DICOM standard application context	1.2.840.10008.3.1.1.1

4.2.2.1 Number of Associations

The service only accepts one incoming connection.

4.2.2.2 Asynchronous Nature

The service does not support asynchronous operations and will not perform asynchronous window negotiation.

4.2.2.3 Implementation Identifying Information

The following fields will be used to identify the implementation:

Field	Value
Implementation class UID	1.3.6.1.4.1.30071.8
Implementation version name	fo-dicom 4.0.4

4.2.3 Association Initiation Policy

The service does not initiate associations.

visi∳∩**rt**

4.2.4 Association Acceptance Policy

The service accepts remote association requests for the following reasons:

- Verify the DICOM protocol communication.
- Store plans from a remote system to the MapRT database.

4.3 Network Interfaces

4.3.1 Physical Network Interface

The MapRT DICOM Listener application works over any physical medium available on the operating system and hardware it is installed on.

4.3.2 Additional Protocols

None.

4.3.3 IPv4 and IPv6 Support

The MapRT DICOM Listener application supports both IPv4 and IPv6.

4.4 Configuration

4.4.1 AE Title/Presentation Address Mapping

4.4.1.1 Local AE Titles.

The application entity title (AET) and TCP/IP port are set in the DICOM server configuration.

4.4.1.2 Remote AE Title/Presentation Address Mapping

Not applicable.

4.4.2 Parameters

The following configuration parameters are set when installing the MapRT DICOM listener:

Name	Default
EntityTitle	MapRT
Port	106



5 MEDIA INTERCHANGE

The MapRT DICOM Listener application does not support Media Interchange.

6 TRANSFORMATION OF DICOM TO CDA

The MapRT DICOM Listener does not output HL7 clinical document architecture reports.



7 SUPPORT OF CHARACTER SETS

The MapRT DICOM Listener application supports all character sets specified by the specific character set attribute (0008,0005) and it is recommended that all DICOM files set this to control how strings are decoded.

If not set, the MapRT DICOM Listener application will default to the UTF7 code page to decode strings. This means that non-ASCII characters, such as accents, can be decoded against a different code page to the one used to encode them, which can result in unexpected behavior such as them being displayed incorrectly.

8 SECURITY

8.1 Security Profiles

No Security Profiles are supported.

8.2 Association Level Security

The MapRT DICOM Listener application will refuse association open requests arriving with an AE title different to that configured.

8.3 Application-Level Security

Local administrative credentials are required to stop, (re)start, or change the configuration of the MapRT DICOM Listener application.

9 ANNEXES

9.1 IOD Contents

9.1.1 Created SOP Instance(s)

No SOP instances are created.

9.1.2 Usage of Attributes from Received IODs

DICOM files stored by the MapRT DICOM Listener can be selected in the following Vision RT product workflows:

9.1.2.1 Clearance Mapping Based on TPS Plan

This process applies to the following SOP classes:

SOP Classes		
Name	UID	
RT Plan	1.2.840.10008.5.1.4.1.1.481.5	

The MapRT DICOM Listener receives RT plans for evaluation from a hospital's treatment planning system (TPS). The patients these refer to will automatically be added to the MapRT database.

It is also possible to manually add a patient by clicking the Add Patient button on the Select Patient screen and selecting a RT plan from the hard disk.

Тад	Name	Description	Notes	
(0008,0005)	Specific Character Set	Character set to expand or replace the basic set.	Optional. Required if an expanded or replacement set is used. Default assumes UTF7.	
(0008,0016)	SOP Class UID	Uniquely identifies the SOP class.	Mandatory and one of: • 1.2.840.10008.5.1.4.1.1.481.5	
(0008,0018)	SOP Instance UID	Uniquely identifies the SOP instance.	Mandatory. Must contain a valid UID (cannot be empty). Duplicate plans (SOP Instance UID already in database) are rejected.	
(0010,0010)	Patient Name	Patient's full name.	Mandatory. Must contain first and last name,	

To load an RT plan the following DICOM tags are used:

Тад	Name	Description	Notes	
			other elements are ignored.	
			If subsequent RT plan with the same Patient ID contains different names than the previous plan imported to database, the plan is rejected.	
	Patient ID	Primary identifier for the Patient.	Mandatory.	
(0010,0020)			Must contain a valid UID (cannot be empty).	
(0010,0030)	Patient Birth Date	Birth date of the Patient.	Mandatory.	
			Cannot be a future date.	
(0010,0040)	Patient Sex	Sex of the patient	Optional.	
			One of 'M' (male), 'F' (female), or 'O' (other). Default is 'O' other.	
	Patient Position	Patient position descriptor relative to the equipment.	Mandatory.	
(0018,5100)			Must be one of HFP, HFS, HFDL, HFDR, FFP, FFS, FFDL or FFDR.	
			Used to convert from DICOM to Vision RT coordinates when later constructing 3D surfaces.	
(300A,0002)	RT Plan Label	User-defined label for treatment plan.	Mandatory.	
			Can be an empty string.	
	Beam Sequence	Sequence of treatment beams for this RT plan.	Mandatory.	
(300A,00B0)			When set, the sequence must describe at least one beam.	

DICOM Conformance Statement

Тад	Name	Description	Notes	
>(300A,00C0)	Beam Number	Identification number of the beam	Mandatory. Must be unique in this plan.	
>(300A,00C2)	Beam Name	User-defined name for beam.	Optional. Default sets an empty string.	
>(300A,00B2)	Treatment Machine Name	User-defined name identifying treatment machine used for treatment delivery.	Optional.	
>(300A,0111)	Control Point Sequence	Sequence of machine configs describing this treatment beam.	Mandatory. When set it must contain at least one control point.	
>>(300A,011E)	Gantry Angle	Gantry angle of the radiation source.	Mandatory in first item of the control point sequence. Optional in subsequent items.	
>>(300A,0122)	Patient Support Angle	Orientation of the patient support coordinate system (couch) in degrees with respect to the fixed reference coordinate system.	Mandatory in the first Item of the control point sequence. Optional in subsequent items.	
>>(300A,011F)	Gantry Rotation Direction	Direction of Gantry Rotation when viewing gantry from isocenter, for segment following Control Point.	Mandatory for the first item of the control point sequence, or if gantry rotation direction changes during beam. One of 'CW' (clockwise), 'CC' (counter-clockwise), 'NONE' (no rotation).	

Tag	Name	Description	Notes		
>>(300A,012C)	Isocenter Position	The 3D co- ordinates (X, Y, Z) of the upper left- hand corner in mm (center of first voxel transmitted).	Mandatory in the first Item of the control point sequence. Optional in subsequent items.		
(300A,0180)	Patient Setup Sequence	Sequence of patient setup data for the current plan.	Mandatory. One or more items shall be included in this Sequence. Used to access patient position (0018,5100) for this beam.		

9.1.3 Attribute Mapping

No attribute mapping is performed.

9.1.4 Coerced/Modified Fields

No attributes are coerced or modified.

9.2 Data Dictionary of Private Attributes

Private Attributes within received files are ignored.

9.3 Coded Terminology and Templates

No Coded Terminology or Templates are supported.

9.4 Grayscale Image Consistency

No additional grayscale adjustments are made when displaying images.

9.5 Standard Extended/Specialized/Private SOP Classes

No Standard Extended, Specialized, or Private SOP Classes are supported.

9.6 Private Transfer Syntaxes

No Private Transfer Syntaxes are supported.

1016-0591-EN (DOC-4579) Ver. 2

Approved By:

Andrzej Wawrzynczyk - Author	
October 17, 2024 12:50 PM BST	286b1b24-a389-424f-9b60-a48ae0834f1c
Amy Mansfield - Regulatory	
October 18, 2024 11:29 AM BST	f65e6f13-a327-4106-b18c-a5e7f154382c
Paul Irwin - Technical	
October 18, 2024 3:39 PM BST	54a7a81b-c2af-41dc-ba0a-4a92db3642a0
Richard Grey - Document Control	
October 22, 2024 3:48 PM BST	37e79ddc-e1ed-4975-a9e8-c984a8f504cf

Version History:

Author	Effective Date	CO#	Ver.	Status
Andrzej Wawrzynczyk	October 22, 2024 3:48 PM BST	Not Available	2	Published
Michelle Kennedy	February 16, 2024 5:46 PM GMT	Not Available	1	Superseded
Nigel Russell	August 8, 2023 4:08 PM BST	<u>CO-93</u>	0	Superseded