



Declaration of Compatibility and Interoperability

Issue Date: 4th July 2025

Vendor Product(s) Validated:

Supported Vendor Product Configurations:

- AlignRT 7.4 , DoseRT module and Patient ID module

Note: DoseRT and Patient ID do not interface to Elekta products and have been assessed for compatibility of the in-room hardware with EMLA only.

Excluded Vendor Product Configurations: None

Elekta Product(s), interfaces and version number(s) as applicable:

Supported Elekta Product Configurations:

- MOSAIQ 2.86 or higher – This includes MOSAIQ in a CITRIX/CLOUD deployment
- The following Elekta Medical Linear Accelerator (EMLA) products
 - Elekta Synergy
 - Elekta Infinity
 - Elekta Harmony
 - Elekta Harmony Pro
 - Elekta Versa HD
 - Elekta Evo

This includes IEC60601-2-1 Edition 3 and Edition 4 compliant EMLA products

Supported Interfaces:

- MOSAIQ PPS – Table moves
- Elekta Response – Beam Hold (EMLA IEC60601-2-1 Edition 3 only)
- Integrated Gating interface - Beam Hold (EMLA IEC60601-2-1 Edition 4 only)

Installation of Integrity 4.1.0.3 is a prerequisite of this Interop declaration

Excluded Elekta Product Configurations:

- Elekta Unity
- Elekta ONE online workflows




Conclusion:

This document certifies that the specified VisionRT Product(s) and the specified Elekta Products have been validated as Interoperable and Compatible.

This declaration applies only to the specified products and versions listed above and shall not be applicable if any changes or updates are made to such specified products and versions when such changes or updates affect the Compatibility and/or Interoperability.

This declaration is valid for 2 years from the Issue Date.

VisionRT Limited

By 

Name NS

Title _____

Date Jul 4, 2025

Elekta Limited

By Christopher Busch

Name CB

Title _____

Date Jul 5, 2025

Elekta Inc

By Anish Patankar

Name AP

Title _____

Date Jul 8, 2025